

STEERING WHEEL HOLDERS



Our steering wheel spinners are designed for reduced hand functions and can be placed either on the right side or on the left side of the steering wheel.




FV197

STEERING WHEEL FORK SPINNER: this is used when the wrist joint mobility is good, but there is a complete lack of mobility in the fingers. All Kivi spinners are easily removable from the steering wheel.



PV197

KNOB SPINNER: It is used in situations where there is only a partial lack of mobility in the wrist and the fingers.



TV197

STEERING WHEEL TRIPOD SPINNER: this is used when there is a complete lack of mobility in both the wrist joint and the fingers.



www.kivi.it | info@kivi.it

facebook.com/Kivi.Fan.page

youtube.com/KiviMobilityFreedom

DEALER